



Innovation for Enhanced Workflow



Harmony iQ[®] 2800
OR Integration System



A Compact Approach to Improving Operative Workflow... Today and Tomorrow

STERIS recognizes the challenges facing today's OR teams:

- Achieving efficient workflow through tools that are easy to use
- Connecting and displaying video from any device
- Keeping focus on the patient to deliver excellent care instead of diagnosing OR equipment issues
- Increasing procedural turnaround and the hospital's bottom line

Through our OR Integration Systems, we focus on supporting patients and hospital teams alike, now and into the future, through easy-to-use and scalable solutions that support the operative workflow.



Connect More: Route 8 Devices to 8 Destinations

The Harmony iQ 2800 Integration System supports the connection of up to 8 devices, routing them to display on up to 8 destinations, as well as capturing and recording in high definition.

Further, it enables HD streaming and conferencing for real-time collaboration through Call and Connect; all from an easy-to-use, centralized control panel.





A User-Focused Touch Panel Designed By Nurses, For Nurses

It's critical to ensure that OR teams can use the technologies within the OR in an easy and efficient manner in order to receive the full benefit of those technologies.

With one, ideally located point of control, the OR staff can operate audio sources, connected computers, in-light cameras, wall cameras and the Free5® camera simply and conveniently. Support efficient and safe procedures, and respond quickly to events that create a need for adjustment to these devices.

More Information. Less Clutter.

Our small, low-profile base unit design connects a wide range of imaging and video devices through a vendor-neutral approach, while maximizing valuable, useable OR space. It easily fits under or on top of a desk.


The unit and cabling are kept off the floor, reducing clutter and trip hazards.



Truly Connected for Better Patient Outcomes

Compare and Contrast Video at a Glance

Experience enhanced visualization that allows OR teams to view video from multiple sources side-by-side on one screen.



Capture Everything. Keep What You Want.

Our **VuCapture® Digital Recording System** enables you to capture surgical procedures in stunning high-definition detail and associate them with patient case* data, then print or save to a USB drive.

Through our **Harmony iQ Perspectives® Image Management System**, easily send captured media to the Electronic Medical Record, and securely review the media from a computer or mobile device. Images and recordings are automatically saved to the correct case to support accurate and thorough case records.*

* Images are not suitable for diagnostic applications.



Happy Patient. Happy Staff.

Through our **SignatureSuite™ OR Integration System** offerings, we can help reduce the patient's anxiety before a procedure by dimming the lights and displaying calming media on monitors, such as a beach scene or mountain landscape, as the patient enters the OR. This feature also helps to reduce stress for the OR team and make set-up run more efficiently.

It's one more tangible way your hospital can demonstrate its focus on the patient. Once the patient is anesthetized, video routing presets exit from the SignatureSuite OR Integration System and prepare for the particular type of procedure's visualization requirements.

Connect and Collaborate in Real-Time.

Connect teams within the OR with clinicians and audiences outside the OR for real-time collaboration and education.

Our **Call and Connect** streaming and conferencing technologies are supported on the Harmony iQ 2800 Integration System, delivering HD audio and video streaming and bi-directional audio and videoconferencing capabilities within the hospital network or beyond.


Ready for 4K? We Are.



4K Planning With Ease and Flexibility

Making your OR 4K compatible can be as easy as upgrading your existing integration equipment.

STERIS' 4K upgrade for the Harmony iQ 2800 Integration System consists of 4K routing hardware that installs with the existing system and the Vividimage® 4K Surgical Displays and Vividimage® W 4K Wall Displays.”



Vendor-Neutral OR Integration: Get the Most Out of Your 4K Investment

What good is a 4K integration system that ties your hands to limited imaging devices and capabilities? Our vendor-neutral architecture empowers you to choose the 4K imaging devices that best suit your needs.

Supporting Your Evolving Technology Needs ... at Your Pace

Our scalable design and fiber-optic architecture support more connections and higher resolutions as your needs grow. Through our future-ready design, you can be assured your OR Integration System will meet your needs today and in the future.

Advanced Visualization Without Compromise

We recognize that 4K OR integration is about more than just seeing ultra high-resolution surgical video. With the Harmony iQ 2800 Integration System, you also retain the robust collaboration and workflow tools you expect in a STERIS integrated OR.





Technology and Support You Can Count On

Skilled STERIS technicians are available by phone and on-site to support real-time identification and resolution of technical issues, as well as perform regular maintenance to keep the system performing at its best. Optional remote diagnostic support lets STERIS perform many support functions without needing to enter the OR.

We'll quickly diagnose and resolve system issues and help new users learn the ropes, keeping your focus on the patient.

STERIS – A History of OR Innovation

STERIS OR Integration Systems support positive clinical outcomes and efficiency, building on generations of surgical solutions that give OR staff the freedom and confidence to focus on what they do best – care for patients.

For more information, contact your STERIS representative or visit www.steris.com.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
800-548-4873
www.steris.com