

Flexible, Efficient and Safe

Harmony® ConnectPoint® Video Connectivity System

Integrate Your OR and Increase the Surgeon's Visibility
for a Fraction of the Cost of a Complete Room Renovation

The Harmony CONNECTPOINT Video Connectivity system offers power and video signal connections from a flexible arm attached to any STERIS surgical lighting system.

It improves workflow and visibility for surgical teams, allowing them to place mobile equipment where they want it; not just where power outlets dictate. Further, the system allows you to:

- Reduce floor trip hazards
- Decrease cleaning time between cases
- Deliver connectivity features of a boom at a lower price

Integrate Your Connections for Safety and Efficiency

Seamlessly Route Video Signals Around the Surgical Field

The Harmony ConnectPoint Video Connectivity System supports all standard video signals (DVI, HD-SDI, S-video and RGBHV), easily delivering surgical video to connected displays or OR integration systems.



The Harmony ConnectPoint Video Connectivity System offers versatile connectivity around the sterile field.



Connect and Position Equipment Where the Surgeon Needs it

The simple, ceiling-mounted arm allows the Harmony ConnectPoint Video Connectivity System to cover a wide 340 degree arc around the sterile field.

Support Quicker, Safer Room Turnaround

The Harmony ConnectPoint Video Connectivity System locates mobile equipment cables off of the floor, removing trip hazards and allowing the surgical staff to perform their tasks more efficiently.



Replace the need for wall outlets and reduce the number of cables routed on the floor, improving safety and efficiency.

For more information, contact your STERIS representative or visit www.steris.com.

 **STERIS**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
800-548-4873
www.steris.com