



MOBILE CART FOR AMSCO® SINGLE COMPARTMENT WARMING CABINET

Durable and functional mobile cart specifically built by Stanley InnerSpace for STERIS's 24" deep single warming cabinet.



Locking casters assure cabinet stays where placed, but allows for easy cleaning around and underneath.

Convenient & Efficient Mobility

Bring adaptability and functionality to the operating or critical care room in the hospital or Ambulatory Surgery Center (ASC). Easily transport the cart from room to room or change positioning within the room for convenient access of warm blankets, IV and Irrigation Fluids, helping to maintain efficiency during procedures.

Sturdy, Secure Construction

Designed specifically for STERIS's 24" Single Chamber Warming Cabinet, this sturdy cart is made with Aluminum Composite Side Panels and 14 GA powder coated steel top and base. Its 1/2" transparent acrylic doors provide easy access and visibility to supplies. The adjustable shelf glides easily into the Stanley InnerSpace patented inner cell system (additional shelves can be purchased separately, part no. XS78SHELF).

Easy-to-grip handles



Technical Specifications

Part No. XS780037WCNC

Width: 34.5"

Depth: 29.5"

Height: 37.5"

Body: Aluminum Composite Panels

Corners: Extruded Aluminum

Base and Top - exterior is made out of 14 GA. Powder coated steel

Casters: 5" Two front locking plate casters with thread guard and zinc plated steel

Doors: 1/2" Thick Transparent Acrylic

Handles: PP plastic

Cell System: ABS plastic

For more information contact your STERIS representative or visit www.steris.com.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com